


Haemostatic Valves/ Angioplasty Pack (Short screw Type)

SUNMED™



Features:

- 1. Maximize the availability of catheters.*
- 2. Minimize the interruption by adapters connected.*
- 3. Suitable for different operations.*

Ordering Information:

Ref (Incomplete)	Description
SM-HV-8SSF-TB	8F, Short screw Type, Sideon Female Luer, Tyvek-plastic Blister
SM-HV-8SSL-TB	8F, Short screw Type, Sideon Tubing and Female Luer, Tyvek-plastic Blister
SM-HV-9SSF-TB	9F, Short screw Type, Sideon Female Luer, Tyvek-plastic Blister
SM-HV-9SSL-TB	9F, Short screw Type, Sideon Tubing and Female Luer, Tyvek-plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 Whatsapp: +1 949 365 6648
Email: sunnymedical_cn@139.com
Website: www.sunnymedical.com www.sunnymedical.com.cn

CE₂₇₉₇
ISO13485
FDA510(k)



Features:

- 1. Ergonomics plunger design provides comfortable palm pushing.*
- 2. Wing grips design suits for various finger sizes of the doctors and nurses.*

Ordering Information:

Ref (Incomplete)	Description
SM-CS-10RPW-TB	10ml, Rotating Male Luer, Palm Pad, Wing Grips, Tyvek-plastic Blister
SM-CS-12RPW-TB	12ml, Rotating Male Luer, Palm Pad, Wing Grips, Tyvek-plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 Whatsapp: +1 949 365 6648
Email: sunnymedical_cn@139.com
Website: www.sunnymedical.com www.sunnymedical.com.cn

CE₂₇₉₇
ISO13485
FDA510(k)
NMPA

Haemostatic Valves/ Angioplasty Pack (Push-screw Type)

SUNMED™



Features:

- 1. Large 10F inner diameter enables two Balloon Catheters or two stent delivery catheters into inner lumen at one time.*
- 2. The use of a hemostatic valve not featured on conventional Y connectors reduces blood loss during procedures.*
- 3. A locking feature of the hemostatic valve has been added, while maintaining the ease-of-use operation.*

Ordering Information:

Ref (Incomplete)	Description
SM-HV-10PSF-TB	10F, Push-screw Type, Sideon Female Luer, Tyvek-plastic Blister
SM-HV-10PSFLCR-TB	10F, Push-screw Type, Sideon Female Luer, Insertion Tool, Red/Copper Torquer, Tyvek-plastic Blister
SM-HV-10PSFLPW-TB	10F, Push-screw Type, Sideon Female Luer, Insertion Tool, White/Plastic Torquer, Tyvek-plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172

Tel: +86 755 8453-2326 Whatsapp: +1 949 365 6648


Email: sunnymedical_cn@139.com

Website: www.sunnymedical.com www.sunnymedical.com.cn

CE₂₇₉₇
ISO13485
FDA510(k)

Disposable Balloon Inflation Devices (J-type)

SUNMED™



Features:

- 1. Easy release mechanism ensures one hand operation.*
- 2. Ergonomic handle design provides comfortable grip.*
- 3. Stable pressure gauge with maximum pressure at 40atm.*

Ordering Information:

Ref (Incomplete)	Description
SM-ID-J2830H-TB	28ml, 30atm/bar, HP-Stopcock, Tyvek-plastic Blister
SM-ID-J2840H-TB	28ml, 40atm/bar, HP-Stopcock, Tyvek-plastic Blister
SM-ID-J2830HP14-TB	28ml, 30atm/bar, HP Stopcock, with P14 Y connector Kit, Tyvek-plastic Blister
SM-ID-J2840HP14-TB	28ml, 40atm/bar, HP Stopcock, with P14 Y connector Kit, Tyvek-plastic Blister
SM-ID-J2830HS14-TB	28ml, 30atm/bar, HP Stopcock, with S14 Y connector Kit, Tyvek-plastic Blister
SM-ID-J2840HS14-TB	28ml, 40atm/bar, HP Stopcock, with S14 Y connector Kit, Tyvek-plastic Blister
SM-ID-J2830HC14-TB	28ml, 30atm/bar, HP Stopcock, with C14 Y connector Kit, Tyvek-plastic Blister
SM-ID-J2840HC14-TB	28ml, 40atm/bar, HP Stopcock, with C14 Y connector Kit, Tyvek-plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172

Tel: +86 755 8453-2326 Whatsapp: +1 949 365 6648

Email: sunnymedical_cn@139.com

Website: www.sunnymedical.com www.sunnymedical.com.cn

CE₂₇₉₇
ISO13485
FDA510(k) Exempt

Haemostatic Valves/ Angioplasty Pack (Y Click)

SUNMED™



Features:

- 1. Three-stage valve creates maneuverability without any bleeding.*
- 2. Auto-closing mechanism ensures the device to move from semi-closed position to closed position automatically when contrast is injected from the side-on lumen.*
- 3. Ergonomic design allows for one hand operation.*

Ordering Information:

Ref (Incomplete)	Description
SM-HV-7CF-TB	7F, Y Click, Sideon Female Luer, Tyvek-plastic Blister
SM-HV-7CP-TB	7F, Y Click, Sideon Tubing and Stopcock, Tyvek-plastic Blister
SM-HV-7CFLCR-TB	7F, Y Click, Sideon Female Luer, Insertion Tool, Red/Copper Torquer, Tyvek-plastic Blister
SM-HV-7CPLCR-TB	7F, Y Click, Sideon Tubing and Stopcock, Insertion Tool, Red/Copper Torquer, Tyvek-plastic Blister
SM-HV-8CFLPW-TB	8F, Y Click, Sideon Female Luer, Insertion Tool, White/Plastic Torquer, Tyvek-plastic Blister
SM-HV-8CPLPW-TB	8F, Y Click, Sideon Tubing and Stopcock, Insertion Tool, White/Plastic Torquer, Tyvek-plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172

Tel: +86 755 8453-2326 Whatsapp: +1 949 365 6648

Email: sunnymedical_cn@139.com

Website: www.sunnymedical.com www.sunnymedical.com.cn

CE₂₇₉₇
ISO 13485

Haemostatic Valves/ Angioplasty Pack (Double Screw Type)

SUNMED™



Features:

- 1. Two separate valves in one adaptor fulfill two guide wires operating independently according to actual needs;*
- 2. Different valves make it easier to distinguish and operate the two guide wires;*
- 3. 2-in-1 design saves procedure time and cost.*

Ordering Information:

Ref (Incomplete)	Description
SM-HV-8WF-TB	8F, Double Screw , Sideon Female Luer, Tyvek-plastic Blister
SM-HV-8WFLCR-TB	8F, Double Screw, Sideon Female Luer, Insertion Tool, Red/Copper 2 Torquers, Tyvek-plastic Blister
SM-HV-8WFLPW-TB	8F, Double Screw, Sideon Female Luer, Insertion Tool, White/Plastic 2 Torquers, Tyvek-plastic Blister
SM-HV-9WF-TB	9F, Double Screw, Sideon Female Luer, Tyvek-plastic Blister
SM-HV-9WFLCR-TB	9F, Double Screw, Sideon Female Luer, Insertion Tool, Red/Copper 2 Torquers, Tyvek-plastic Blister
SM-HV-9WFLPW-TB	9F, Double Screw, Sideon Female Luer, Insertion Tool, White/Plastic 2 Torquers, Tyvek-plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 Whatsapp: +1 949 365 6648
Email: sunnymedical_cn@139.com
Website: www.sunnymedical.com www.sunnymedical.com.cn

CE₂₇₉₇
ISO13485
FDA510(k)

Haemostatic Valves/ Angioplasty Pack (Push-pull Type)

SUNMED™



Features:

- 1. Two different sideons of Y adaptor suit for various customers.*
- 2. Push-pull Y connector makes it easier for doctor to operate precisely.*

Ordering Information:

Ref (Incomplete)	Description
SM-HV-7PF-TB	7F, Push-pull, Sideon Female Luer, Tyvek-plastic Blister
SM-HV-8PP-TB	8F, Push-pull, Sideon Tubing and Stopcock, Tyvek-plastic Blister
SM-HV-7PPLCR-TB	7F, Push-pull, Sideon Tubing and Stopcock, Insertion Tool, Red/Copper Torquer, Tyvek-plastic Blister
SM-HV-8PFLCR-TB	8F, Push-pull, Sideon Female Luer, Insertion Tool, Red/Copper Torquer, Tyvek-plastic Blister
SM-HV-7PFLPW-TB	7F, Push-pull, Sideon Female Luer, Insertion Tool, White/Plastic Torquer, Tyvek-plastic Blister
SM-HV-8PPLPW-TB	8F, Push-pull, Sideon Tubing and Stopcock, Insertion Tool, White/Plastic Torquer, Tyvek-plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 Whatsapp: +1 949 365 6648
Email: sunnymedical_cn@139.com
Website: www.sunnymedical.com www.sunnymedical.com.cn

CE2797
ISO13485
FDA510(k)
NMPA

Haemostatic Valves/ Angioplasty Pack (Screw Type)

SUNMED™



Features:

- 1. Two different types of Y adaptor suit for various customers.*
- 2. Large closing knob makes twisting and valve closing very easy.*

Ordering Information:

Ref (Incomplete)	Description
SM-HV-8SF-TB	8F, Screw Type, Sideon Female Luer, Tyvek-plastic Blister
SM-HV-8SP-TB	8F, Screw Type, Sideon Tubing and Stopcock, Tyvek-plastic Blister
SM-HV-8SFLCR-TB	8F, Screw Type, Sideon Female Luer, Insertion Tool, Red/Copper Torquer, Tyvek-plastic Blister
SM-HV-8SPLCR-TB	8F, Screw Type, Sideon Tubing and Stopcock, Insertion Tool, Red/Copper Torquer, Tyvek-plastic Blister
SM-HV-8SFLPW-TB	8F, Screw Type, Sideon Female Luer, Insertion Tool, White/Plastic Torquer, Tyvek-plastic Blister
SM-HV-8SPLPW-TB	8F, Screw Type, Sideon Tubing and Stopcock, Insertion Tool, White/Plastic Torquer, Tyvek-plastic Blister
SM-HV-9SF-TB	9F, Screw Type, Sideon Female Luer, Tyvek-plastic Blister
SM-HV-9SP-TB	9F, Screw Type, Sideon Tubing and Stopcock, Tyvek-plastic Blister
SM-HV-9SFLCR-TB	9F, Screw Type, Sideon Female Luer, Insertion Tool, Red/Copper Torquer, Tyvek-plastic Blister
SM-HV-9SPLCR-TB	9F, Screw Type, Sideon Tubing and Stopcock, Insertion Tool, Red/Copper Torquer, Tyvek-plastic Blister
SM-HV-9SFLPW-TB	9F, Screw Type, Sideon Female Luer, Insertion Tool, White/Plastic Torquer, Tyvek-plastic Blister
SM-HV-9SPLPW-TB	9F, Screw Type, Sideon Tubing and Stopcock, Insertion Tool, White/Plastic Torquer, Tyvek-plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.


1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172

Tel: +86 755 8453-2326 Whatsapp: +1 949 365 6648

Email: sunnymedical_cn@139.com

Website: www.sunnymedical.com www.sunnymedical.com.cn

CE2797
ISO13485
FDA510(k)



Features:

- 1. Luer lock fits both ordinary syringe and control syringe.*
- 2. Tape with super adhesive guarantees air balloon having sufficient stanching pressure.*

Ordering Information:

Ref (Incomplete)	Description
SM-AP-L-MP	Large size 30cmx26cm, Medical paper-plastic pouch


Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 Whatsapp: +1 949 365 6648
Email: sunnymedical_cn@139.com
Website: www.sunnymedical.com www.sunnymedical.com.cn

CE₂₇₉₇
ISO13485
FDA510(k)
NMPA

Disposable Balloon Inflation Devices (M-type)

SUNMED™



Features:

- 1. Spiral design ensures precise pressure.*
- 2. Improved release handle makes instant and easy inflation and deflation.*


Ordering Information:

Ref (Incomplete)	Description
SM-ID-M2030-TB	20ml, 30atm/bar, Tyvek-plastic Blister
SM-ID-M2030H-TB	20ml, 30atm/bar, HP Stopcock, Tyvek-plastic Blister
SM-ID-M2030S01-TB	20ml, 30atm/bar, with S01 Y connector Kit, Tyvek-plastic Blister
SM-ID-M2030P01-TB	20ml, 30atm/bar, with P01 Y connector Kit, Tyvek-plastic Blister
SM-ID-M2030C01-TB	20ml, 30atm/bar, with C01 Y connector Kit, Tyvek-plastic Blister
SM-ID-M2030W01-TB	20ml, 30atm/bar, with W01 W connector Kit, Tyvek-plastic Blister
SM-ID-M2030PS01-TB	20ml, 30atm/bar, with PS01 Y connector Kit, Tyvek-plastic Blister
SM-ID-M2030HS01-TB	20ml, 30atm/bar, HP Stopcock, with S01 Y connector Kit, Tyvek-plastic Blister
SM-ID-M2030HP01-TB	20ml, 30atm/bar, HP Stopcock, with P01 Y connector Kit, Tyvek-plastic Blister
SM-ID-M2030HC01-TB	20ml, 30atm/bar, HP Stopcock, with C01 Y connector Kit, Tyvek-plastic Blister
SM-ID-M2030HW01-TB	20ml, 30atm/bar, HP Stopcock, with W01 W connector Kit, Tyvek-plastic Blister
SM-ID-M2030HPS01-TB	20ml, 30atm/bar, HP Stopcock, with PS01 Y connector Kit, Tyvek-plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 Whatsapp: +1 949 365 6648
Email: sunnymedical_cn@139.com
Website: www.sunnymedical.com www.sunnymedical.com.cn

CE₂₇₉₇
ISO13485
FDA510(k) Exempt
NMPA



Features:

1. Blue dot leads to precise positioning of TR-closure Band.
2. Three different sizes are applicable to almost all patients. Large size, Medium size and Small size.
3. Pressure in air balloons is easily adjusted by specially designed syringe tip adaptor.

Ordering Information:

Ref (Incomplete)	Description
SM-TB-L-TP	Large Size without Syringe, Tyvek-plastic Pouch
SM-TB-S-TP	Small Size without Syringe, Tyvek-plastic Pouch
SM-TB-LS-TP	Large Size with 20ml Syringe, Tyvek-plastic Pouch
SM-TB-SS-TP	Small Size with 20ml Syringe, Tyvek-plastic Pouch
SM-TB-LLS-TP	Large Size with Syringe, Luer Connector, Tyvek-plastic Pouch
SM-TB-MLS-TP	Medium Size with Syringe, Luer Connector, Tyvek-plastic Pouch
SM-TB-SLS-TP	Small Size with Syringe, Luer Connector, Tyvek-plastic Pouch

Sunny Medical Device (Shenzhen) Co., Ltd.


1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172

Tel: +86 755 8453-2326 Whatsapp: +1 949 365 6648

Email: sunnymedical_cn@139.com

Website: www.sunnymedical.com www.sunnymedical.com.cn

CE²⁷⁹⁷
ISO13485
NMPA



Features:

- 1. Specially designed 2ml deflation lock mechanism facilitates accurate air discharge.*
- 2. Central balloon design stops bleeding quickly.*
- 3. Bandage on the belt holds the air syringe in place for the doctor to operate.*
- 4. Three different sizes are applicable to almost all patients. Large size, Medium size and Small size.*

Ordering Information:

Ref (Incomplete)	Description
SM-TB-CLLS-TP	C Type, Large Size with Syringe, Luer Connector, Tyvek-plastic Pouch
SM-TB-CMLS-TP	C Type, Medium Size with Syringe, Luer Connector , Tyvek-plastic Pouch
SM-TB-CSLS-TP	C Type, Small Size with Syringe, Luer Connector , Tyvek-plastic Pouch


Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172

Tel: +86 755 8453-2326 Whatsapp: +1 949 365 6648

Email: sunnymedical_cn@139.com

Website: www.sunnymedical.com www.sunnymedical.com.cn



Features:

- 1. Rotating male luer provides convenient and fast connection.*
- 2. Lubricated rotating handles enhance handle rotatability.*

Ordering Information:

Ref (Incomplete)	Description
SM-MF-25WF2B-TB	250psi, White Handles, Off Handles, 2 Ports, Block Body, Tyvek-plastic Blister
SM-MF-25WO2B-TB	250psi, White Handles, On Handles, 2 Ports, Block Body, Tyvek-plastic Blister
SM-MF-25BO3H-TB	250psi, Blue Handles, On Handles, 3 Ports, Half Body, Tyvek-plastic Blister
SM-MF-25BF3BP-TB	250psi, Blue Handles, Off Handles, 3 Ports, Block Body, Protectors, Tyvek-plastic Blister
SM-MF-50BO3B-TB	500psi, Blue Handles, On Handles, 3 Ports, Block Body, Tyvek-plastic Blister
SM-MF-50BF3B-TB	500psi, Blue Handles, Off Handles, 3 Ports, Block Body, Tyvek-plastic Blister
SM-MF-50WO3B-TB	500psi, White Handles, On Handles, 3 Ports, Block Body, Tyvek-plastic Blister
SM-MF-50WF3B-TB	500psi, White Handles, Off Handles, 3 Ports, Block Body, Tyvek-plastic Blister
SM-MF-50BF3BP-TB	500psi, Blue Handles, Off Handles, 3 Ports, Block Body, Protectors, Tyvek-plastic Blister

Note: All manifolds are right-handed.

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 Whatsapp: +1 949 365 6648
Email: sunnymedical_cn@139.com
Website: www.sunnymedical.com www.sunnymedical.com.cn

CE₂₇₉₇
ISO13485
FDA510(k) Exempt
NMPA



100psi



400psi
Swivel Nut



500psi



1200psi

Features:

- 1. Stopcock with rotating male luer can meet high pressure application requirement.*


Ordering Information:

Ref (Incomplete)	Description
SM-SC-04BFSF-MP	400psi, Blue, Off handle, Swivel Nut, Female luer, Medical paper plastic Pouch
SM-SC-05BFRF-TB	500psi, Blue Handle, Off Handle, Rotating Male Luer, Female Luer, Tyvek-Plastic Blister
SM-SC-05WORF-TB	500psi, White Handle, On Handle, Rotating Male Luer, Female Luer, Tyvek-Plastic Blister
SM-SC-12BFRF-TB	1200psi, Blue Handle, Off Handle, Rotating Male Luer, Female Luer, Tyvek-Plastic Blister
SM-SC-12WORF-TB	1200psi, White Handle, On Handle, Rotating Male Luer, Female Luer, Tyvek-Plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 Whatsapp: +1 949 365 6648
Email: sunnymedical_cn@139.com
Website: www.sunnymedical.com www.sunnymedical.com.cn

CE₂₇₉₇
ISO13485
FDA510(k) Exempt
NMPA



Features:

- 1. Saving costs by all-in-one package.*
- 2. Saving time with connected and ready to use kit.*

Ordering Information:

Ref (Incomplete)	Description
SM-MK-FI2C1P1-TP	Manifold with Off and White Handle, Block Body, Control Syringe, 2 Infusion Sets, 1 Pressure Line, Tyvek-plastic Pouch
SM-MK-FI2C1P1-PP	Manifold with Off and White Handle, Block Body, Control Syringe, 2 Infusion Sets, 1 Pressure Line, Plastic-plastic Pouch

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 Whatsapp: +1 949 365 6648
Email: sunnymedical_cn@139.com
Website: www.sunnymedical.com www.sunnymedical.com.cn



Features:

- 1. Transparent cylinder enhances visual effect of gauge mark and contrast level.*
- 2. Special designed Silicone piston head and lubrication technology lead plunger easy to move in cylinder and generate high pressure resistance.*
- 3. Thumb and finger rings fullfill precise dose control.*

Ordering Information:

Ref (Incomplete)	Description
SM-CS-10RRR-TB	10ml, Rotating Male Luer, Thumb Ring, Ring Grips, Tyvek-plastic Blister
SM-CS-12RRR-TB	12ml, Rotating Male Luer, Thumb Ring, Ring Grips, Tyvek-plastic Blister
SM-CS-12RRRL-TB	12ml, Rotating Male Luer, Thumb Ring, Ring Grips, Locking Plunger, Tyvek-plastic Blister
SM-CS-10FRR-TB	10ml, Fixed Male Luer, Thumb Ring, Ring Grips, Tyvek-plastic Blister
SM-CS-12FRR-TB	12ml, Fixed Male Luer, Thumb Ring, Ring Grips, Tyvek-plastic Blister
SM-CS-12FRRL-TB	12ml, Fixed Male Luer, Thumb Ring, Ring Grips, Locking Plunger, Tyvek-plastic Blister


Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 Whatsapp: +1 949 365 6648
Email: sunnymedical_cn@139.com
Website: www.sunnymedical.com www.sunnymedical.com.cn

CE₂₇₉₇
ISO13485
FDA510(k)
NMPA

Pressure Lines

SUNMED™



Coiled Pressure Line



Features:

1. 400psi/500psi/900psi/1200psi pressure line by PVC Tubing.
2. Low cost, high performance.
3. 1200psi Pressure Line by Braided PU Tubing.
4. Flexible tubing and rotating male luer provide easy and fast connection.


Ordering Information:

Ref (Incomplete)	Description
SM-PL-04PC120MFP-MP	400psi, PVC, Coiled, 120cm, Male to Female, Protectors, Medical paper-plastic Pouch
SM-PL-04PC120MF-MP	400psi, PVC, Coiled, 120cm, Male to Female, Medical paper-plastic Pouch
SM-PL-05P120MF-MP	500psi, PVC, 120cm, Male to Female, Medical paper-plastic Pouch
SM-PL-12P120MF-MP	1200psi, PVC, 120cm, Male to Female, Medical paper-plastic Pouch
SM-PL-12P120RF-MP	1200psi, PVC, 120cm, Rotating Male to Female, Medical paper-plastic Pouch
SM-PL-12B120RF-TP	1200psi, Braided PU, 120cm, Rotating Male to Female, Tyvek-plastic Pouch

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 Whatsapp: +1 949 365 6648
Email: sunnymedical_cn@139.com
Website: www.sunnymedical.com www.sunnymedical.com.cn

CE₂₇₉₇
ISO13485
FDA510(k)
NMPA



Features:

1. *sharp tip that has low resistance increases accessibility and reduce tissue trauma.*
2. *Short tip has better controllability.*
3. *Tapered hub inner diameter facilitates insertion of the guide wire.*
4. *Needles of customized sizes are also available.*


Ordering Information:

Ref (Incomplete)	Description
SM-PN-18U07LS-MP	18G, Ultra-Thin Wall, 7cm, Large Hub, Short Tip, Medical Paper-plastic Pouch
SM-PN-18U07SM-MP	18G, Ultra-Thin Wall, 7cm, Small Hub, Medium Tip, Medical Paper-plastic Pouch
SM-PN-18U07BM-MP	18G, Ultra-Thin Wall, 7cm, B Hub, Medium Tip, Medical Paper-plastic Pouch
SM-PN-18T07LM-TB	18G, Thin Wall, 7cm, Large Hub, Medium Tip, Tyvek-plastic Blister
SM-PN-18T07SM-TB	18G, Thin Wall, 7cm, Small Hub, Medium Tip, Tyvek-plastic Blister
SM-PN-18T07BS-TB	18G, Thin Wall, 7cm, B Hub, Short Tip, Tyvek-plastic Blister
SM-PN-21T04LM-MP	21G, Thin Wall, 4cm, Large Hub, Medium Tip, Medical Paper-plastic Pouch
SM-PN-21U04LM-MP	21G, Ultra-Thin Wall, 4cm, Large Hub, Medium Tip, Medical Paper-plastic Pouch
SM-PN-21T04LM-TB	21G, Thin Wall, 4cm, Large Hub, Medium Tip, Tyvek-plastic Blister
SM-PN-21U04LM-TB	21G, Ultra-Thin Wall, 4cm, Large Hub, Medium Tip, Tyvek-plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
 Shenzhen, Guangdong, China 518172
 Tel: +86 755 8453-2326 Whatsapp: +1 949 365 6648
 Email: sunnymedical_cn@139.com
 Website: www.sunnymedical.com www.sunnymedical.com.cn

CE2797
 ISO13485
 FDA510(k) Exempt
 NMPA



Features:

1. *3cm flexible tip makes tip easy to deflect.*
2. *Elastic and shape-resistance distal core and lubricious coating provide great accessibility and excellent steerability.*
3. *Supportive proximal shaft increases pushability.*

Ordering Information:

Ref (Incomplete)	Description
SM-GW-3S32150P-TP	3cm flexible S-tip, 0.032"x150cm, PTFE Coated, Tyvek-plastic Pouch
SM-GW-3J32150P-TP	3cm flexible J-tip, 0.032"x150cm, PTFE Coated, Tyvek-plastic Pouch
SM-GW-3S32260P-TP	3cm flexible S-tip, 0.032"x260cm, PTFE Coated, Tyvek-plastic Pouch
SM-GW-3J32260P-TP	3cm flexible J-tip, 0.032"x260cm, PTFE Coated, Tyvek-plastic Pouch
SM-GW-3S35150P-TP	3cm flexible S-tip, 0.035"x150cm, PTFE Coated, Tyvek-plastic Pouch
SM-GW-3J35150P-TP	3cm flexible J-tip, 0.035"x150cm, PTFE Coated, Tyvek-plastic Pouch
SM-GW-3J35180P-TP	3cm flexible J-tip, 0.035"x180cm, PTFE Coated, Tyvek-plastic Pouch
SM-GW-3S35260P-TP	3cm flexible S-tip, 0.035"x260cm, PTFE Coated, Tyvek-plastic Pouch
SM-GW-3J35260P-TP	3cm flexible J-tip, 0.035"x260cm, PTFE Coated, Tyvek-plastic Pouch
SM-GW-3S38150P-TP	3cm flexible S-tip, 0.038"x150cm, PTFE Coated, Tyvek-plastic Pouch
SM-GW-3J38150P-TP	3cm flexible J-tip, 0.038"x150cm, PTFE Coated, Tyvek-plastic Pouch
SM-GW-3S38260PN-TP	3cm flexible S-tip, 0.038"x260cm, PTFE Coated, with N dispenser Tyvek-plastic Pouch
SM-GW-3J38260PN-TP	3cm flexible J-tip, 0.038"x260cm, PTFE Coated, with N dispenser Tyvek-plastic Pouch
SM-GW-3J35150H-TP	3cm flexible J-tip, 0.035"x150cm, Hydrophilic Coated, with N dispenser Tyvek-plastic Pouch
SM-GW-3J35180H-TP	3cm flexible J-tip, 0.035"x180cm, Hydrophilic Coated, with N dispenser Tyvek-plastic Pouch
SM-GW-3J35260H-TP	3cm flexible J-tip, 0.035"x260cm, Hydrophilic Coated, with N dispenser Tyvek-plastic Pouch

Note: every J-tip wire has a 3mm radius "J".


Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
 Shenzhen, Guangdong, China 518172
 Tel: +86 755 8453-2326 Whatsapp: +1 949 365 6648
 Email: sunnymedical_cn@139.com
 Website: www.sunnymedical.com www.sunnymedical.com.cn

CE₂₇₉₇
 ISO13485
 FDA510(k)
 NMPA



Guide wire dispenser straightener



Guide wire dispenser adaptor

Guide wire without coating



Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172

Tel: +86 755 8453-2326 Whatsapp: +1 949 365 6648

Email: sunnymedical_cn@139.com

Website: www.sunnymedical.com www.sunnymedical.com.cn

CE₂₇₉₇
ISO13485
FDA510(k)



Ordering Information:

Ref (Incomplete)	Description
SM-HC-T-MP	Transparent, Medical Paper-plastic Pouch

disposable towel/sheet clips



Ordering Information:

Ref (Incomplete)	Description
SM-TC-W1-MP	White, 1 in 1 Pouch, Medical Paper-plastic Pouch
SM-TC-W2-MP	White, 2 in 1 Pouch, Medical Paper-plastic Pouch
SM-TC-W5-MP	White, 5 in 1 Pouch, Medical Paper-plastic Pouch

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172

Tel: +86 755 8453-2326 Whatsapp: +1 949 365 6648

Email: sunnymedical_cn@139.com

Website: www.sunnymedical.com www.sunnymedical.com.cn

CE₂₇₉₇
ISO13485
NMPA



Rotating male luer



Fixed male luer



Fixed female luer



Small hub



Large hub



B hub




Large hub with wing

Clipless Dispenser



Stopper



Protector



Sunny Medical Device (Shenzhen) Co., Ltd.


1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172

Tel: +86 755 8453-2326 Whatsapp: +1 949 365 6648

Email: sunnymedical_cn@139.com

Website: www.sunnymedical.com www.sunnymedical.com.cn

CE₂₇₉₇
ISO13485



Features:

- 1. Superior visualization of fluid;*
- 2. Consistent and accurate monitoring;*
- 3. Easy assembly, reduce setup time;*
- 4. High quality chip from global leading sensor supplier,
Measurement Specialities.*

Ordering Information:

Ref (Incomplete)	Description
SM-DT-UT1-TB	1 Transducer with Utah Type of Cable Connector (TCC), Tyvek-plastic Blister
SM-DT-BB1-TB	1 Transducer with B Braun TCC, Tyvek-plastic Blister
SM-DT-BD1-TB	1 Transducer with BD TCC, Tyvek-plastic Blister
SM-DT-ED1-TB	1 Transducer with Edward TCC, Tyvek-plastic Blister
SM-DT-AB1-TB	1 Transducer with Abbott TCC, Tyvek-plastic Blister
SM-DT-UT1AP2S1-TB	1 Transducer with Utah TCC, 1 Administration Set, 2 Pressure Lines, One Stopcock, Tyvek-plastic Blister
SM-DT-BB1AP2S1-TB	1 Transducer with B Braun TCC, 1 Administration Set, 2 Pressure Lines, One Stopcock, Tyvek-plastic Blister
SM-DT-BD1AP2S1-TB	1 Transducer with BD TCC, 1 Administration Set, 2 Pressure Lines, One Stopcock, Tyvek-plastic Blister
SM-DT-ED1AP2S1-TB	1 Transducer with Edward TCC, 1 Administration Set, 2 Pressure Lines, One Stopcock, Tyvek-plastic Blister
SM-DT-AB1AP2S1-TB	1 Transducer with Abbott TCC, 1 Administration Set, 2 Pressure Lines, One Stopcock, Tyvek-plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.


1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172

Tel: +86 755 8453-2326 Whatsapp: +1 949 365 6648

Email: sunnymedical_cn@139.com

Website: www.sunnymedical.com www.sunnymedical.com.cn

Plate/Organizer



Clamp/Mounting Bracket



Male-Female Luer



Adapter Cable



Utah Interface



Abbott Interface



Ordering Information:

Ref (Incomplete)

Description

OEM-DT-P-NBP	Plate for Transducer
OEM-DT-C-NBP	Clamp for Transducer
OEM-DT-AC-UT2MR	IBP Adapter Cable Utah Pressure Transducer to Mindray monitor.
OEM-DT-AC-UT2HP	IBP Adapter Cable Utah Pressure Transducer to Philips/HP monitor.
OEM-DT-AC-AB2MR	IBP Adapter Cable Abbott Pressure Transducer to Mindray monitor.
OEM-DT-AC-AB2HP	IBP Adapter Cable Abbott Pressure Transducer to Philips/HP monitor.

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District, Shenzhen, Guangdong, China 518172

Tel: +86 755 8453-2326 Whatsapp: +1 949 365 6648


Email: sunnymedical_cn@139.com

Website: www.sunnymedical.com www.sunnymedical.com.cn

Disposable Balloon Inflation Devices

Digital display pressure gauge (A-type)

SUNMED™



Features:

- 1. Ergonomic design facilitates doctor's operation.*
- 2. Trigger release mechanism provides easy and quick deflation.*
- 3. Inclination of manometer provides easy and precise view of pressure points.*
- 4. Digital display pressure gauge for precise control of operational pressure.*

Ordering Information:

Ref (Incomplete)	Description
SM-ID-AD2530H-TB	25ml, 30atm/bar(Digital display), HP-Stopcock, Tyvek-plastic Blister
SM-ID-AD2530HP14-TB	25ml, 30atm/bar(Digital display), HP Stopcock, with P14 Y connector Kit, Tyvek-plastic Blister
SM-ID-AD2530HS14-TB	25ml, 30atm/bar(Digital display), HP Stopcock, with S14 Y connector Kit, Tyvek-plastic Blister
SM-ID-AD2530HC14-TB	25ml, 30atm/bar(Digital display), HP Stopcock, with C14 Y connector Kit, Tyvek-plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172

Tel: +86 755 8453-2326 Whatsapp: +1 949 365 6648


Email: sunnymedical_cn@139.com

Website: www.sunnymedical.com www.sunnymedical.com.cn

CE₂₇₉₇
ISO13485

Disposable Balloon Inflation Devices (S-type) 30ATM(40ATM)

SUNMED™



Features:

- 1. Square shape handle applies to different sizes of hands.*
- 2. Chinese knots on handle prevent slippery on holding.*
- 3. Soft cover on deflation release button reduces deflation vibration.*

Ordering Information:

Ref (Incomplete)	Description
SM-ID-S3030H-TB	S-type 30ml, 30atm/bar, HP Stopcock, Tyvek-plastic Blister
SM-ID-S3040H-TB	S-type 30ml, 40atm/bar, HP Stopcock, Tyvek-plastic Blister
SM-ID-S3030HP14-TB	S-type 30ml, 30atm/bar, HP Stopcock, with P14 Y Connector Kit Tyvek-plastic Blister
SM-ID-S3030HS14-TB	S-type 30ml, 30atm/bar, HP Stopcock, with S14 Y Connector Kit, Tyvek-plastic Blister
SM-ID-S3030HC14-TB	S-type 30ml, 30atm/bar, HP Stopcock, with C14 Y Connector Kit, Tyvek-plastic Blister
SM-ID-S3040HP14-TB	S-type 30ml, 40atm/bar, HP Stopcock, with P14 Y Connector Kit Tyvek-plastic Blister
SM-ID-S3040HS14-TB	S-type 30ml, 40atm/bar, HP Stopcock, with S14 Y Connector Kit, Tyvek-plastic Blister
SM-ID-S3040HC14-TB	S-type 30ml, 40atm/bar, HP Stopcock, with C14 Y Connector Kit, Tyvek-plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172

Tel: +86 755 8453-2326 Whatsapp: +1 949 365 6648


Email: sunnymedical_cn@139.com

Website: www.sunnymedical.com www.sunnymedical.com.cn

CE₂₇₉₇
ISO13485
FDA510(k)

Disposable Balloon Inflation Devices (A-type)

SUNMED™



Features:

- 1. Ergonomic design facilitates doctor's operation.*
- 2. Trigger release mechanism provides easy and quick deflation.*
- 3. Inclination of manometer provides easy and precise view of pressure points.*

Ordering Information:

Ref (Incomplete)	Description
SM-ID-A2530H-TB	25ml, 30atm/bar, HP-Stopcock, Tyvek-plastic Blister
SM-ID-A2530HP14-TB	25ml, 30atm/bar, HP Stopcock, with P14 Y connector Kit, Tyvek-plastic Blister
SM-ID-A2530HS14-TB	25ml, 30atm/bar, HP Stopcock, with S14 Y connector Kit, Tyvek-plastic Blister
SM-ID-A2530HC14-TB	25ml, 30atm/bar, HP Stopcock, with C14 Y connector Kit, Tyvek-plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172
Tel: +86 755 8453-2326 Whatsapp: +1 949 365 6648
Email: sunnymedical_cn@139.com
Website: www.sunnymedical.com www.sunnymedical.com.cn

CE₂₇₉₇
ISO13485
NMPA

Disposable Balloon Inflation Devices (F-type)

SUNMED™



Features:

- 1. Accurate inflation and safe deflation;*
- 2. Ergonomic design facilitates doctor's operation;*
- 3. Specially designed automatic cam locking makes it easy-handling.*

Ordering Information:

Ref (Incomplete)	Description
SM-ID-F3030H-TB	F-type, 30ml, 30atm/bar, HP Stopcock, Tyvek-plastic Blister
SM-ID-F3030HP14-TB	F-type, 30ml, 30atm/bar, HP Stopcock, with P14 Y connector Kit, Tyvek-plastic Blister
SM-ID-F3030HS14-TB	F-type, 30ml, 30atm/bar, HP Stopcock, with S14 Y connector Kit, Tyvek-plastic Blister
SM-ID-F3030HC14-TB	F-type, 30ml, 30atm/bar, HP Stopcock, with C14 Y connector Kit, Tyvek-plastic Blister

Sunny Medical Device (Shenzhen) Co., Ltd.


1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172

Tel: +86 755 8453-2326 Whatsapp: +1 949 365 6648

Email: sunnymedical_cn@139.com

Website: www.sunnymedical.com www.sunnymedical.com.cn

CE₂₇₉₇
ISO13485



- **Ideal for CTO**
- **Semi-Compliance**
- **6F Guide Kissing**

MediSunik's semi-compliance pre-dilatation catheter which offers numerous features for enhanced crossability, flexibility and extraordinary rewrap to facilitate the physician and procedure. The balloon sizes 1.0 and 1.25mm are designed for small and narrow lesions and for the treatment of CTOs. It is no boundary to cross and no limitation, GO FURTHER AND GO FASTER.

- *Manicured and balanced 0.016" tip*
- *Enhanced balloon flexibility and crossability*
- *5F guiding catheter compatibility and 6F guiding catheter compatibility for kissing balloon technique*
- *Innovative shaft materials with unsurpassed kink resistance*

Balloon Specification

Balloon Diameter(mm)	1.0, 1.25, 1.5, 1.75, 2.0, 2.25, 2.5, 2.75, 3.0, 3.25, 3.5, 3.75, 4.0
Balloon Length(mm)	5, 8, 10, 12, 15, 20, 25, 30
Nominal Pressure	6atm
Rated Burst Pressure	14atm
Average Burst Pressure	20atm
Compliance	Semi-compliance
Balloon Pleats	2/3
Balloon Profile	0.021"(ø1.0mm); 0.030"(ø3.0mm)
Tip Length(mm)	1.5/2.0/2.5
Distal Shaft OD	1.9F
Proximal Shaft OD	2.36F/2.55F/2.7F
Guiding Catheter Compatibility	5F Compatible/6F compatible for kissing
Working Length(cm)	140
Coating	Slider™ Hydro from tip to adjacent guide wire port (external surface); Eel™ coating for guide wire lumen


Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172

Tel: +86 755 8453-2326 Whatsapp: +1 949 365 6648

Email: sunnymedical_cn@139.com

Website: www.sunnymedical.com www.sunnymedical.com.cn



• High Pressure • Non-Compliance • 6F Guide Kissing

MediSunik's truly non-compliance high pressure post-dilatation catheter which offers numerous features for enhanced robust & non-compliance balloon, extraordinary minor balloon axial growth while inflation and 6F guiding catheter compatible. The balloon sizes 2.0-5.0mm and length 8-18mm. The average balloon burst pressure is 30atm. It is satisfied and unlimited, the STRONGER the BETTER and SAFER.

- Manicured and balanced 0.016" tip
- Nominal @ 12atm and Rated Burst Pressure @ 22atm
- Minimal Balloon Compliance, 0.55% from nominal to Rated Burst Pressure(RBP)
- Axial Balloon Growth, 3% from Nominal to Rated Burst Pressure
- Average burst pressure 30atm
- 5F guiding catheter compatibility and 6F guiding catheter compatibility for kissing balloon technique
- Innovative shaft materials with unsurpassed kink resistance

Balloon Specification

Balloon Diameter(mm)	2.0, 2.25, 2.5, 2.75, 3.0, 3.25, 3.5, 3.75, 4.0, 4.5, 5.0
Balloon Length(mm)	8, 10, 12, 15, 18
Nominal Pressure	12atm
Rated Burst Pressure	22atm(∅2.0-4.0mm) 20atm(∅4.5-5.0mm)
Average Burst Pressure	30atm
Compliance	Non-compliance, 0.55%
Balloon Pleats	3
Tip Length(mm)	2.0/2.5
Distal Shaft OD	2.55/2.6F
Proximal Shaft OD	2.0F
Guiding Catheter Compatibility	5F Compatible/6F compatible for kissing
Working Length(cm)	140
Coating	Slider™ Hydro from tip to adjacent guide wire port (external surface); Eel™ coating for guide wire lumen

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172

Tel: +86 755 8453-2326 Whatsapp: +1 949 365 6648

Email: sunnymedical_cn@139.com

Website: www.sunnymedical.com www.sunnymedical.com.cn



Ordering Information:

Ref (Incomplete)	Description
SM-FS-001	Transparent, 330mm*220mm, 1 in 1 Pouch, Non-sterile, Bulk Package

Sunny Medical Device (Shenzhen) Co., Ltd.

1/F and 401, Zhongtianxin Building B, Longgang District,
Shenzhen, Guangdong, China 518172

Tel: +86 755 8453-2326 Whatsapp: +1 949 365 6648

Email: sunnymedical_cn@139.com

Website: www.sunnymedical.com www.sunnymedical.com.cn

ISO13485

FDA510(k) Exempt